

UROMED CYSTOBAG TK 1000 S 30/S 60/S 90

geschlossenes Urindrainagesystem zur Tag- und Nachtversorgung

Closed Urine Drainage System for Day and Night Care

REF 4870/4871/4872

1000 ml

Schlauchlänge 30 cm/60 cm/90 cm
Vliesbeschichtung

1000 ml

tube length 30 cm/60 cm/90 cm
fleece cover

Der UROMED CYSTOBAG TK 1000 für die Tag- und Nachtversorgung vereint die Vorteile eines geschlossenen Urindrainagesystems mit den Vorzügen eines Beinbeutels und der damit verbundenen Mobilität des Patienten. Ein häufiges Diskonnektieren verschiedener Systeme am Kathetertrichter entfällt. Zur Befestigung am Ober- bzw. Unterschenkel empfehlen wir unser UROMED Klettbandset REF 4896.


UROMED CYSTOBAG TK 1000 for day and night care combines the advantages of a closed urine drainage system with the features of an urine leg bag and the mobility of the patient involved. A frequent disconnecting of different systems at the catheter funnel will be stopped. For the fixation at the femur or lower leg we recommend our UROMED hook-and-loop-strap set REF 4896.

Steriler Urinbeutel zum Auffangen von Urin nach Ableitung durch einen Katheter

- mit latexfreien Handschuhen
- Stufenkonnektor für zuverlässige Anschlüsse
- nadelfreie Urinprobeentnahmestelle
- Tropfkammer mit Rücklaufventil
- Belüftung an Tropfkammer und Beutel
- Fixationskordel zur universellen Befestigung
- übersichtliche Graduierung
- T-TAP-Ablasshahn Halterung
- T-TAP-Ablasshahn einhändig und kontaminationsfrei bedienbar
- Vliesbeschichtung auf der Beutelnrückseite

Sterile urine bag for collection of urine drained with a catheter

- with gloves free of latex
- step connector for reliable links
- needle-free spot for urine specimen
- drip chamber with non-return valve
- air vent on drip chamber and bag
- strip for all-purpose fixation
- clearly graduated
- dockhouse
- T-TAP outlet one-handed usable without contamination
- fleece cover on the reverse side of the bag



REF	Schlauchlänge	Volumen	PZN
REF	Tube length	Volume	CPN
4870	30 cm	1000 ml	00689326
4871	60 cm	1000 ml	00689332
4872	90 cm	1000 ml	00689349

Hilfsmittelkatalog: 15.25.07.1003

Verpackung: steril
Originalkarton: 20 Stück

Packaging: sterile
Original packing: 20 pieces